

Comparison of Full-Scan, Data-Dependent, and Data-Independent Acquisition Modes in Liquid Chromatography-Mass Spectrometry Based Untargeted Metabolomics

Jian Guo, and Tao Huan

Anal. Chem., **Just Accepted Manuscript** • DOI: 10.1021/acs.analchem.9b05135 • Publication Date (Web): 13 May 2020

Downloaded from pubs.acs.org on May 14, 2020

Just Accepted

“Just Accepted” manuscripts have been peer-reviewed and accepted for publication. They are posted online prior to technical editing, formatting for publication and author proofing. The American Chemical Society provides “Just Accepted” as a service to the research community to expedite the dissemination of scientific material as soon as possible after acceptance. “Just Accepted” manuscripts appear in full in PDF format accompanied by an HTML abstract. “Just Accepted” manuscripts have been fully peer reviewed, but should not be considered the official version of record. They are citable by the Digital Object Identifier (DOI®). “Just Accepted” is an optional service offered to authors. Therefore, the “Just Accepted” Web site may not include all articles that will be published in the journal. After a manuscript is technically edited and formatted, it will be removed from the “Just Accepted” Web site and published as an ASAP article. Note that technical editing may introduce minor changes to the manuscript text and/or graphics which could affect content, and all legal disclaimers and ethical guidelines that apply to the journal pertain. ACS cannot be held responsible for errors or consequences arising from the use of information contained in these “Just Accepted” manuscripts.

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

**Comparison of Full-Scan, Data-Dependent, and Data-Independent Acquisition Modes in
Liquid Chromatography-Mass Spectrometry Based Untargeted Metabolomics**

Jian Guo¹, Tao Huan^{1*}

¹ Department of Chemistry, Faculty of Science, University of British Columbia, Vancouver
Campus, 2036 Main Mall, Vancouver, V6T 1Z1, BC, Canada

* Corresponding author: Dr. Tao Huan, thuan@chem.ubc.ca

Abstract

Full-scan, data-dependent acquisition (DDA), and data-independent acquisition (DIA) are the three common data acquisition modes in high resolution mass spectrometry-based untargeted metabolomics. It is an important yet underrated research topic on which acquisition mode is more suitable for a given untargeted metabolomics application. In this work, we compared the three data acquisition techniques using a standard mixture of 134 endogenous metabolites and a human urine sample. Both hydrophilic interaction and reversed-phase liquid chromatographic separation along with positive and negative ionization modes were tested. Both the standard mixture and urine samples generated consistent results. Full-scan mode is able to capture the largest number of metabolic features, followed by DIA and DDA (53.7% and 64.8% respective features fewer on average in urine than full-scan). Comparing the MS² spectra in DIA and DDA, spectra quality is higher in DDA with average dot product score 83.1% higher than DIA in Urine(H), and the number of MS² spectra (spectra quantity) is larger in DIA (on average 97.8% more than DDA in urine). Moreover, a comparison of relative standard deviation distribution between modes shows consistency in the quantitative precision, with the exception of DDA showing a minor disadvantage (on average 19.8% and 26.8% fewer features in urine with RSD < 5% than full-scan and DIA). In terms of data preprocessing convenience, full-scan and DDA data can be processed by well-established software. In contrast, several bioinformatic issues remain to be addressed in processing DIA data and the development of more effective computational programs is highly demanded.

31 Introduction

32 Untargeted metabolomics is emerging as a powerful and attractive technology that has been widely
33 used in many research areas, including, but not limited to, the discovery of disease biomarkers,
34 study of metabolic pathways, monitoring the totality of environmental exposures, and searching
35 for bioactive chemicals that modulate certain disease phenotypes, such as demyelinating disease,
36 diabetes, cancer, etc.¹⁻⁵. Liquid chromatography coupled to high resolution mass spectrometry
37 (LC-HRMS) is currently the most prominent analytical platform for untargeted metabolomics,
38 owing to its high throughput, sensitivity and specificity. In LC-HRMS based untargeted
39 metabolomics, LC separates metabolites by their chemical properties to reduce sample complexity,
40 and retention time information of metabolites can be used to assist their identification. Following
41 separation, MS is used to detect metabolic signal intensities that are proportional to metabolite
42 concentrations for quantitative assessments. There are three MS data acquisition modes to record
43 metabolic signals: full-scan, data-dependent acquisition (DDA), and data-independent acquisition
44 (DIA). Conventionally, LC-HRMS based untargeted metabolomics is carried out by first
45 performing LC-full scan MS on each individual sample to obtain accurate mass-to-charge ratio
46 (m/z) and relative abundance for all the metabolic features. The acquired metabolic information is
47 then processed using data preprocessing software (*e.g.*, XCMS, MZmine)^{6,7} to automatically
48 extract high-confidence metabolic features. Further multivariate and univariate statistical analyses
49 extract significantly altered metabolic features that are pertinent to the given biological question⁸.
50 From there, LC-MS² (also called LC-MS/MS or LC-tandem MS) analysis is performed on the
51 pooled sample to collect fragmentation spectra of metabolic features of interest in order to confirm
52 their chemical identities. This conventional approach has been well established but it separates
53 metabolite quantification from structural determination, thus requires extra instrumental analysis
54 time and sample volumes.

55

56 To facilitate a more time-efficient approach of data acquisition in untargeted metabolomics, an

57 accelerated and autonomous data acquisition mode, DDA is implemented. In DDA mode, the MS

58 instrument performs MS full-scan immediately followed by MS² analysis on a list of precursor

59 ions selected from the full-scan spectrum⁹. In this setting, both quantitative (obtained from the MS

60 full-scan) and structural (obtained the MS² spectra) information can be acquired in the same sample

61 analysis, allowing for simultaneous data processing and metabolite identification. Since the

62 selection of precursor ions for MS² analysis is intensity dependent, a potential limitation of DDA

63 mode is that metabolic features of interest with low MS abundance may never be selected for

64 fragmentation and will have no MS² data associated with them. In addition, since the MS

65 instrument allocates a large portion of its acquisition time for MS² spectra generation, the signal

66 intensity for MS¹ features is reduced. This causes a fundamental problem for the detection and

67 quantification of low abundance metabolites in DDA mode.

68

69 To address the limitation of DDA in acquiring MS² spectra, recent implementation of DIA in

70 metabolomic studies can theoretically generate MS² spectra for all precursor ions, therefore

71 capable of detecting and identifying more metabolites that are at lower concentrations^{10,11}. In DIA

72 mode, the MS instrument cycles through the precursor ion m/z range with a large precursor ion

73 mass width to fragment more than one precursor ion simultaneously. Commonly used DIA

74 methods include all-ion fragmentation (AIF) (*e.g.*, MS^{All}. MS^E)^{12,13}, where all precursors are

75 fragmented, and sequential window acquisition of all theoretical fragment-ion spectra (SWATH)¹⁴,

76 where a medium pass window, such as 20 or 25 Da, is used. Theoretically, DIA enables continuous

77 and unbiased acquisition of MS² information for all metabolites. Although LC-HRMS in DIA

mode is an attractive metabolomic profiling strategy, the link between precursors and their fragment ions are dissociated due to the complexity of the resulting MS² spectra. Thus a great informatics challenge to deconvolute the raw MS² spectra for subsequent metabolite identification is posed¹⁵.

Given the unique features and limitations of each data acquisition mode, it is rational to believe that the choice of data acquisition mode may affect the quality of metabolomics results. It is also important to note that a common practice in targeted studies is to maximize the detection sensitivity via optimizing sample preparation and LC-MS conditions. However, it is much harder in untargeted analysis to do the abovementioned optimization due to the diverse chemical properties and concentrations of all metabolites. Therefore, it is more practical to select the best data acquisition mode for better analytical performance in untargeted metabolomics. Unfortunately, there is no direct comparison of the three data acquisition modes yet, and the literature of choosing an appropriate acquisition mode for given untargeted metabolomics projects is inadequate. At the current stage, most researchers are more inclined to use the acquisition mode that has already been established in their lab without considering which acquisition method might work better. To address this knowledge gap, we made a comparative study of full-scan, DDA, and DIA AIF modes in LC-HRMS based untargeted metabolomics. The capability of different data acquisition modes was tested using metabolite standard mixture and human urine samples. Complementary LC-MS configurations were assessed, including hydrophilic interaction liquid chromatography (HILIC) and reversed-phase liquid chromatography (RP), as well as electrospray ionization (ESI) in positive and negative modes. Metabolomic data collected in the three data acquisition modes were evaluated from five complementary aspects: (1) metabolic coverage; (2) quantitative precision; (3)

1
2
3 101 MS² spectral quality; (4) MS² spectral coverage, the percentage of metabolic features that have
4
5 102 MS² spectra out of the total metabolic features, and (5) convenience. Our results demonstrate that
6
7 103 full-scan provides the best feature detection sensitivity and quantitative precision. In comparison,
8
9 104 DDA provides high quality MS² spectra and DIA provides better MS² spectral coverage.
10
11
12

13 105
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

106 **Experimental Section**

107 **Chemicals and solvents.** LC-MS grade water (H₂O), acetonitrile (ACN), Methanol (MeOH) and
108 other solvents were purchased from Fisher Scientific (Hampton, New Hampshire). All the
109 metabolite standards and other chemicals were purchased from Sigma-Aldrich (St. Louis,
110 Missouri).

111
112 **Standard solution preparation.** The 134 metabolite standards were first prepared individually as
113 stock solutions, each in 1 mL ACN:H₂O (1:1, v:v). To prepare the original standard mixture, all
114 the metabolites were pooled together to reach a final concentration of 8.3 µg/mL. To prepare
115 diluted standard mixture, 10 µL original solution was further diluted 100 times to 8.3 ng/mL. The
116 original standard mixture is indicated as “StdMix(H)” and the diluted standard mixture is indicated
117 as “StdMix(L)”. The list of the 134 metabolites is available in **Table S-1**.

118
119 **Urine sample preparation.** One urine sample was collected from a healthy male volunteer. 50
120 µL urine was mixed with 150 µL LC-MS grade MeOH in a 1.5 mL Eppendorf vial. After vortexing,
121 the solution was spun down at 17,530 g, and the clear urine supernatant was aliquoted to a clean
122 1.5 mL Eppendorf vial. The clear solution was evaporated using speed-vac, and the dried urine
123 sample was reconstituted in 50 µL ACN:H₂O (1:1, v:v) for LC-HRMS analysis. A diluted urine
124 sample was prepared by a 100 time dilution of a 10 µL aliquot from the reconstituted solution. The
125 reconstituted urine solution before dilution is indicated as “Urine(H)” and the diluted solution is
126 indicated as “Urine(L)”.

LC-HRMS analysis. The LC-HRMS analysis was performed on a Bruker Impact II™ Ultra-High Resolution Qq-Time-Of-Flight Mass Spectrometer (UHR-QqTOF-MS)) coupled with an Agilent 1290 Infinity™ II Ultra High-Performance Liquid Chromatography (UHPLC) system. The resolution of this instrument is up to 50,000 (measured at full width at half maximum (FWHM) and at m/z 1222 (positive mode) and m/z 1334 (negative mode) using TuneMix solution and factory default methods) and the mass accuracy is below 5 ppm. 2 μ L sample solutions of each triplicate and solvent blank were injected in sequence onto a Waters ACQUITY UPLC BEH C18 column (130Å, 1.7 μ m, 1.0 mm \times 100 mm) in RP mode, and a Millipore ZIC-pHILIC column (200Å, 5 μ m particle size, 2.1 \times 150 mm) in HILIC mode. 2 μ L of 150 mM Sodium Formate (NaFA) solution was injected for internal calibration.

We tested the data acquisition modes in RPLC separation coupled with electrospray ionization (ESI) in positive mode and HILIC separation coupled with ESI in negative mode. The combination of RP(+) and HILIC(-) is a well-established method in metabolomics to achieve a global-scale metabolome profiling^{16,17}. In RP analysis, mobile phase A was H₂O in 0.1% Formic Acid (FA) and mobile phase B was ACN in 0.1% FA. RP gradient was set as follows: 0 min, 95% A; 8 min, 75% A; 14 min, 30% A; 20 min, 5% A; 23 min, 5% A. Post-gradient equilibrium time was 7 min at 95% A. Flow rate was 0.15 mL/min. In HILIC analysis, mobile phase A was 95% H₂O and 5% ACN with 10 mM Ammonium Acetate at 9.8 pH, and mobile phase B was 5% H₂O and 95% ACN. HILIC gradient was set as follows: 0 min, 5% A; 20 min, 80% A; 25 min, 95% A. Post-gradient equilibrium time was 19 min at 5% A. Flow rate was 0.15 mL/min during the gradient, and 0.2 mL/min during the post-gradient equilibrium. Q-TOF MS was performed in full-scan, DDA, and DIA AIF modes. The detailed MS parameter settings can be found in **Text S-1**.

151

152 **Data processing.** MS-DIAL was used for metabolomic data processing. The computer had an

153 Intel i9-9900k CPU @ 3.60 GHz with 8 cores and 32 GB memory; Windows 10; 64-bit operation

154 system; and 10 processing threads. The parameters were optimized for processing full-scan, DDA

155 and DIA data^{18,19}. The detailed parameters can be found in **Text S-2**. After metabolic feature

156 extraction, statistical analysis was performed to calculate the relative standard deviations of

157 metabolic signals in analytical replicates (**Text S-3**). All the raw data has also been uploaded on

158 Metabolights, a public-accessible metabolomics data repository²⁰. Researchers can access the raw

159 data via its job ID: MTBLS1572.

160

161 **Metabolite annotation.** For the standard mixture, the metabolites were annotated by matching

162 their accurate mass and retention time to a tandem MS library file created in house using the same

163 instrument and pure chemical standards, thereby containing the accurate mass and retention time

164 for each metabolite in the standard mixture. For urine samples, the metabolites were annotated by

165 MS² matching against the publicly available MS spectral database provided by MS-DIAL (MSMS-

166 Public-Pos-VS11 and MSMS-Public-Neg-VS11). The library contains accurate mass and tandem

167 MS information for annotation. The positive mode library includes the tandem MS information for

168 29,269 metabolites and the negative mode library includes the tandem MS information for 17,810

169 metabolites. The tandem MS information was collected from Massbank, ReSpec, GNPS, Fiehn

170 lab, CASMI2016, MetaboBase, and RIKEN PlaSMA. The instruments used to generate the MS²

171 spectra are Orbitrap and Q-TOF MS. In addition, a MS² quality score was calculated for the MS²

172 spectral matching by comparing the number of common product ions between the evaluated

173 spectrum and the reference spectrum from MS-DIAL. The score can be represented by a well-

174 established value of dot product. The detailed explanation of dot product calculation can be found
175 in **Text S-4**.

176

177 Results

178 The schematic workflow for the comprehensive assessment of full-scan, DDA, and DIA AIF is
179 illustrated in **Figure 1**. A standard mixture consisting of 134 endogenous metabolites and a healthy
180 human urine sample (both original concentration and 100-fold dilution) were tested in this study.
181 LC-MS analysis was performed in both RP(+) and HILIC(-), the two most commonly used LC-
182 MS modes for profiling hydrophobic and hydrophilic metabolites, respectively^{16,17}. After data
183 collection, the three data acquisition modes were compared in terms of metabolic feature numbers,
184 quantitative precision, MS² spectral quality and quantity, as well as convenience.

185
186 **Analysis of standard mixture.** All the metabolites included in the standard mixture and their
187 detectability in different concentrations and acquisition modes are listed in **Table S-1, 2 and 3**,
188 respectively. The retention time and peak intensity of each identified metabolite in each condition
189 are listed in **Table S-4 and 5**. After manual checking of retention time, accurate mass and MS²
190 spectra of the metabolite standards, 66, 63, and 66 metabolites were detected in full-scan, DDA,
191 and DIA modes in RP analysis, respectively. In HILIC analysis, 125, 120, and 125 metabolites
192 were detected in full-scan, DDA, and DIA modes (**Table 1**), respectively. The metabolites not
193 extracted by MS-DIAL but recovered by manual checking are also listed in **Table S-7**. From these
194 results, we noticed that the number of confirmed metabolites were consistently slightly higher in
195 full-scan and DIA than in DDA in both RP(+) and HILIC(-) modes.

196
197 **Metabolite detection at lower sample concentrations.** The ability to detect low abundance
198 metabolites is sought after in untargeted metabolomics. To assess the performance of different data
199 acquisition modes in detecting low abundance metabolites, we repeated the analysis of metabolite
200 standard mixture using its 100-fold diluted solution (**Table 1**). Applying the same LC-HRMS

analysis to StdMix(L), we confirmed 44, 23, and 33 metabolites in full-scan, DDA and DIA modes in RP(+) analysis and 57, 43, and 46 metabolites in full-scan, DDA, and DIA modes in HILIC(-) analysis. As expected, with 100-fold dilution, the number of detected metabolites was significantly reduced compared with undiluted original solution (StdMix(H)). Similar to the results obtained using the StdMix(H), full-scan performed the best in terms of the number of detected metabolites. However, compared with StdMix(H) results, DDA had the greatest decrease in the number of identified metabolites (**Table 1**), indicating that DDA is less sensitive in detecting low abundance metabolites than full-scan and DIA. The detailed information about the detected metabolites in StdMix(H) and StdMix(L) are listed in **Table S-1**.

MS² spectra quality. In addition to the quantity of detected metabolites, we also compared the quality of MS² spectra generated in DDA and DIA modes using the same metabolite standards. The collected MS² spectra were compared against the reference MS² spectra archived in the MS-DIAL MS² spectra library. Dot product was calculated in MS-DIAL to assess the similarity between experimental and reference MS² spectra. The detailed matching scores are presented in **Table S-6**. The averaged matching score was 72.8 for DDA results and 46.8 for DIA results. Further t-test analysis showed that there is a significant difference in the MS² spectra quality (p-value: 6.4×10^{-17}) between DDA and DIA results, and DDA generates higher quality MS² spectra than DIA.

Analysis of urine samples. In addition to the targeted analysis using metabolite standards, we further performed an untargeted analysis on human urine samples in full-scan, DDA, and DIA modes. After the urine metabolomic data was generated, MS-DIAL was chosen as the data

processing software for metabolic feature extraction, alignment, and preliminary metabolite identification. Data processing parameters were manually optimized to achieve the best performance of metabolic feature extraction. Our results showed that the average processing times for three replicate sample files plus one solvent blank file were: 110 seconds for full-scan, 40 seconds for DDA, and 7,620 seconds (127 minutes) for DIA. Although the processing of DIA data can be faster if the MS¹ and MS² abundance cut-offs were set higher in MS-DIAL, higher cut-offs may result in the loss of some low abundance metabolic features.

To evaluate the analytical difference between the three data acquisition modes, we first counted the total numbers of detected, MS¹ annotated and MS² matched metabolic features (**Figure 2**). To calculate statistical significance among the results of the three data acquisition modes, we also processed the metabolomic raw data individually and therefore have three replicates for each acquisition mode to perform a t-test (**Figure S-1**). Our definitions for these different features are as follows: detected features were the direct output of peak picking software after filtering out low abundance noise, background signals from blank, and different adduct forms; MS¹ annotated features matched *m/z* values with theoretical *m/z* values of at least one metabolite in the MS-DIAL metabolic library; MS² matched features have experimental MS² spectra matched with reference MS² spectra with a minimum of 40% matching fragment ions and 15% similarity of dot product. As shown in **Figure 2 and Figure S-1**, in both RP(+) and HILIC(-) analyses, full-scan and DIA modes were able to cover more metabolic features (statistically significant, *p*-value < 0.05) than DDA mode regardless of the urine sample concentration. Confidently identified metabolites with dot product > 0.8 from Urine(H) samples in both RP(+) and HILIC(-) modes are listed in **Table S-8**.

247

248 **Comparison of Quantitative Precision.** Besides the number of detected metabolic features, we

249 also investigated the quantitative precision using experimental triplicates of the same urine sample.

250 After feature alignment in MS-DIAL, relative standard deviation (RSD) of the metabolic feature

251 intensities were calculated. As shown in **Figure 3**, although all the three datasets have > 80%

252 metabolic features of RSD < 20%, more metabolic features have excellent precision (RSD < 5%)

253 in full-scan and DIA. This trend is consistent in Urine(H) analyzed in both RP(+) and HILIC(-).

254 The percentages of metabolic features with RSD < 20% are similar between acquisition modes.

255 However, in Urine(L) RP(+) mode, Urine(H) RP(+) mode, Urine(L) HILIC(-) mode, and Urine(H)

256 HILIC(-) mode, the respective percentages of metabolic features with RSD < 5% are 26%, 42%,

257 37% and 52% in DDA; 37%, 53%, 39%, and 67% in full-scan; and 40%, 72%, 49%, and 57% in

258 DIA. DDA consistently has the smallest percentage of metabolic features with RSD < 5%,

259 independent of sample type, concentration, ionization polarity or chromatographic method. It is

260 also noticed in Urine(L) analysis that metabolic features have broader RSD distribution, and more

261 metabolic features have larger RSD values in comparison to Urine(H) in both RP(+) and HILIC(-)

262 analyses.

263

264 **Comparison of urine MS² spectra.** Using human urine samples, we also compared the quantity

265 and quality of MS² spectra generated in DDA and DIA modes. The number of metabolic features

266 with MS² spectra were 981, 668, 62, and 151 in Urine(H)-RP(+), Urine(H)-HILIC(-), Urine(L)-

267 RP(+), and Urine(L)-HILIC(-) analyses using DDA mode, respectively. On average, 43% of

268 metabolic features have corresponding MS² spectra in DDA mode. In comparison, the number of

269 metabolic features with MS² spectra were 1381, 1549, 167, and 225 in Urine(H)-RP(+), Urine(H)-

HILIC(-), Urine(L)-RP(+), and Urine(L)-HILIC(-) analyses using DIA mode, respectively. Since all the metabolic MS² spectra were deconvoluted mathematically from raw MS² data, 100% of the metabolic features have corresponding MS² spectra in DIA analysis. Furthermore, by overlapping common metabolic features in DDA and DIA using their retention time and accurate *m/z*, Venn diagrams were prepared to show the shared metabolic features with MS² spectra (**Figure 4**). As shown, DIA was able to generate more MS² matched features than DDA. However, we noticed that DDA also generated a few unique MS² features which were not found in the MS² matched features in DIA mode. We manually investigated these unique features and discovered that these are features with low quality peak shapes or low abundance that triggered MS² events in DDA mode, and was not picked up in DIA mode. The details can be found in the supplementary information (**Text S-5**).

Besides the comparison of MS² spectral numbers, we also evaluated the quality of MS² spectra generated from DDA and DIA analyses. The metabolic features generated in Urine(H) by both RP(+) and HILIC(-) analyses were used to compare the MS² spectra quality between DDA and DIA. Same as in the metabolite standards analysis, the MS² spectra quality of each metabolic feature was assessed by their dot products against reference MS² spectra in MS-DIAL. The dot product score distributions of all common metabolic features detected by DDA and DIA are presented in **Figure 5**. Due to the diverse metabolic concentration, the MS² spectral quality distribution was similar to those of metabolite standards, and several metabolites had extremely low matching dot product score. A comparison of the dot product score distributions for DDA and DIA results shows that MS² spectra generated in DDA have better scores and thus better spectra quality.

293

294 Discussion

295 In this work, a comprehensive assessment of full-scan, DDA, and DIA AIF modes in untargeted
296 metabolomics was performed. We tested the modes using standard mixture and human urine
297 samples in original and diluted concentrations and RP(+) and HILIC(-) configurations. Our results
298 demonstrate that each data acquisition mode has its unique features as well as limitations that can
299 significantly affect the metabolomic data quality. **Figure 6** summarizes the performance ranking
300 of full scan, DDA, and DIA in metabolic feature detection, quantitative precision, MS² spectral
301 quality and quantity, and convenience. The rankings have three levels represented by three grey
302 pentagon-shape lines in the radar plot (**Figure 6**). The five aspects for comparison are indicated
303 by the five arrow dimensions towards five directions. The further a mode is placed from the center
304 of the plot, the better performance it has in the corresponding dimension. Overall, full-scan is
305 highlighted by its capability to detect the largest number of metabolic features and provide the
306 highest quantification precision. In comparison, DDA is convenient to use, as it streamlines
307 metabolic profiling and metabolite identification, but the number of metabolic features and
308 quantification precision are not as good. Furthermore, between DDA and DIA, the two modes
309 capable of generating MS² spectra, DDA provides better quality and DIA provides better quantity.
310 The purpose of this work is to draw the attention of the scientific community to recognize that the
311 choice of data acquisition mode can lead to different metabolomic results. The comprehensive
312 comparison also provides guidance on how to choose the optimal data acquisition mode in order
313 to fully exploit the analytical power of different data acquisition modes in untargeted
314 metabolomics.

1
2
3 316 In this work, the comparison of the three data acquisition modes was carried out in optimized scan
4
5 317 mode parameters using a Bruker Impact II™ UHR-QqTOF-MS. QqTOF-MS is one of the most
6
7 318 commonly used HRMS for untargeted metabolomics. Its fast scan speed and short accumulation
8
9 319 time allows it to routinely perform full-scan, DDA, and DIA (both AIF and SWATH-type)
10
11 320 analyses with high acquisition frequency and good mass resolution. There are three major reasons
12
13 321 for the choice of AIF (or MS^{all}) as the DIA method instead of SWATH or SWATH-like workflow
14
15 322 for this comparison study. First, AIF is the original DIA method which has been commonly used
16
17 323 and widely available for all MS instruments. In comparison, SWATH is only available in several
18
19 324 major MS instruments. There remains to be a number of MS instruments that are not yet capable
20
21 325 of performing SWATH. Therefore, we decided to compare the more accessible DIA method,
22
23 326 which is AIF. More importantly, there are many parameters in SWATH mode that can influence
24
25 327 the data acquisition results. These parameters need to be optimized depending on the instrument.
26
27 328 However, the optimization has been reported in only a few studies²¹, and is a research topic that
28
29 329 needs further research. Last but not least, a previous study comparing DDA, SWATH and AIF for
30
31 330 drug metabolomics has already shown that SWATH should be better than AIF in terms of MS²
32
33 331 spectral quality due to smaller isolation window and thus less complex MS² spectra¹⁰. Although
34
35 332 we did not include SWATH in our study, it is clear that the implementation and optimization of
36
37 333 SWATH in more MS instruments should be addressed in future research²².
38
39
40
41
42
43
44
45
46

47 335 Our results indicate that full-scan mode is able to detect the largest number of metabolic features
48
49 336 compared to DDA and DIA. The large difference in feature numbers can be explained by the time
50
51 337 the instrument spends on MS¹ data acquisition. Detailed acquisition time information can be found
52
53 338 in **Table S-2**. Since in full-scan mode, 100% of the instrument time is spent on MS¹ data
54
55
56
57
58
59
60

339 acquisition, more low abundance metabolic features can be captured, as evidenced in both the
340 analysis of metabolite standards and urine samples. Thus, full-scan mode could be more suitable
341 to achieve largest metabolic coverage. In comparison, DDA and DIA is unable to achieve that
342 level of comprehensiveness because part of the instrument time is used for acquiring MS²
343 information. Specifically, for the methods used in this study, DDA and DIA respectively spent
344 only approximately 5% and 50% of instrument time on MS¹, comparing to 100% in full-scan mode
345 (**Table S-9**). Therefore, the numbers of detected metabolic features in DDA and DIA modes are
346 much lower than full-scan mode. A detailed illustration of how LC-MS chromatogram is
347 constructed in each acquisition mode is presented in **Figure S-2**. Manual inspection of the DDA
348 data, which contains the fewest number of metabolic features, shows that many of the real
349 metabolic features were not picked up by data processing algorithms because of their poor
350 chromatographic peak shapes and low MS signal intensities. These real metabolic features could
351 potentially be retrieved by our recent effort of integrated feature extraction strategy, which
352 recognizes these low abundance peaks by their valid MS² information, rather than their
353 chromatographic peak shapes²³.

354
355 Generating precise quantitative information is important for comparative metabolomics, and
356 therefore technical variations should be minimized during sample preparation and data collection²⁴.
357 In different data acquisition modes, the quantitative precision is predominantly influenced by peak
358 abundance and peak shape. Our results showed that full-scan is able to have the best quantitative
359 precision, followed by DIA and DDA which generates the least quantitative precision. Although
360 more than 80% of the metabolic features generated in these three data acquisitions have RSD
361 below 20%, which is commonly accepted for LC-MS based metabolomics²⁵, we noticed that full-

scan mode is able to generate the best quantitative precision. We also manually checked some metabolic features with large RSD values. These features are noticed to have low MS scan numbers, and the peak shapes of these metabolic features are not smooth enough to achieve the most appropriate peak heights. Therefore, the peak shapes of these metabolic features in DDA are not as feasible as full-scan and DIA for peak height determination, and lead to more features with larger RSD. However, the disadvantage of DDA in terms of quantitative precision is not dramatic, so the ranking of DDA is placed between the first and second level in our summary radar plot (**Figure 6**). Comparatively, as a result of the high quantitative precision and presence of structural information DIA mode has been commonly used for targeted quantification in proteomics²⁶⁻²⁸ and now starts to be used for targeted quantitative analysis in metabolomics. Recently, there has been an informatics tool developed to automate the quantification process²⁹.

MS² spectra generated in LC-HRMS analysis provides valuable structural information for metabolite annotation via either *de novo* interpretation or matching against reference MS² spectra. Besides comparing MS¹ results, we also investigated the quality and quantity of MS² spectra generated in DDA and DIA modes. Since metabolic structure elucidation heavily relies on the quality of its MS² spectrum, higher spectral quality is favored as it improves the matching accuracy, and leads to a lower chance of false positive identification. Our results showed that DDA mode consistently generates MS² spectra with better quality than DIA mode. This result is similar to previous study where DDA and DIA modes were compared for the identification of drug metabolites¹⁰. It is not difficult to understand the result as in DDA experiment, a narrow precursor ion isolation window (approx. 1 Da) is used to select precursor ions for MS². Thus, by isolation of the pure precursor, a “clearer” MS² spectra can be achieved. In contrast, in DIA experiment, MS²

spectra were generated via mathematical deconvolution of a mixture of MS² fragments. For co-eluting metabolites, it is often difficult to allocate fragment signals to the right metabolic precursor ions, and better deconvolution algorithms are highly warranted to generate cleaner/higher quality MS² spectra in DIA. Although DDA mode outcompetes DIA in terms of MS² spectra quality, DIA is able to generate MS² spectra for all the metabolic features. To address the limited number of MS² spectra generation in DDA analysis, several efforts have been made to improve MS² coverage in DDA analysis³⁰⁻³².

It is important to note that convenience is a major consideration when carrying out untargeted metabolomics. In terms of metabolic feature extraction, many well-established bioinformatic programs for full-scan and DDA data are available to automatically extract metabolic features and align them across multiple samples³³⁻³⁵. In particular, DDA is the most convenient method among all three alternatives, as it can be fully automated from feature extraction to metabolite identification and shortens the full-scan-based workflow from days to hours⁹. More importantly, recent development of spectral similarity analysis (e.g., GNPS) also makes DDA the method of choice to generate large number of MS² spectra in each sample for unknown metabolite annotation³⁶. In our study, the processing time for DIA-based data is almost 100 times longer than DDA and full-scan-based data under the same parameter settings. In addition, lots of false positive metabolic features were observed in the low signal-to-noise (S/N) range of the DIA results. Although many DIA deconvolution algorithms have been developed in the past few years³⁷⁻⁴², we believe that the further development of software for processing DIA data will promote the increased use of DIA for untargeted metabolomics.

Among the limitations of the present study is the limited choice of LC-MS combinations. We focused on the most commonly used combinations in LC-MS based metabolomics (e.g., RP(+) and HILIC(-)), which is also common practice in other MS method comparison works^{10,43-46}. Researchers who want to check their specific LC-MS methods can follow our workflow to perform their comparison. Furthermore, the parameters for LC-MS analysis and data processing have been considered optimized, yet the choice is fairly subjective. We performed parameter optimization via (1) checking recommended settings provided by the data processing program; (2) testing different parameter combinations and manually checking which combination gave the best results; (3) comparing our parameter settings with literature for further optimization. Notably, method optimization via strict Design of Experiments approach of a wide parameter space might generate results with stronger scientific rationale⁴⁷. Nevertheless, our results show that screening through a certain space of the key parameters does not cause significant alterations to the overall trend (**Figure S-3**), suggesting that our method optimization approach is acceptable. Finally, mobile phase pH, column temperature and other factors can affect LC chromatographic separation and thus impact data acquisition performance. These parameters comprise an interesting research topic not addressed in and is beyond the scope of this study.

In summary, this work presents a systematic approach to comparing data acquisition modes in untargeted metabolomics. Our results show that all the three data acquisition modes have distinctive analytical performance in metabolic coverage, quantitative precision, MS² spectra quality and quantity, and convenience. Although the underlying difference is simply caused by how instrument time is allocated for MS¹ and MS², the consequence is substantial, leading to

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

430 different analytical performance. Therefore, the selection of an appropriate data acquisition
431 method should be carefully weighed based on the purpose of the study.
432

Supporting Information

The Supporting Information is available free of charge.

Text S-1: Detailed MS parameter settings for full-scan, DDA and DIA modes. Text S-2: Detailed parameters for processing full-scan, DDA and DIA data in MS-DIA. Text S-3: Statistical analysis used in this comparative study. Text S-4: Details of dot product calculation. Text S-5: The reasons of DDA having unique MS² matched features undetected in DIA. Table S-1: List of metabolites in the standard mixture. Table S-2: Metabolites detected in StdMix(H). Table S-3: Metabolites detected in StdMix(L). Table S-4: Retention time and peak intensity information of metabolites detected in StdMix(H). Table S-5: Retention time and peak intensity information of metabolites detected in StdMix(L). Table S-6: Dot product of the common metabolic feature MS² between RP(+) DDA and DIA methods from StdMix(H). Table S-7: Comparison of the statistical results generated using the MS/MS precursor ions and that using manual examination of the chromatogram in the raw data. Table S-8: Table of identified metabolites in Urine(H) samples in both RP(+) and HILIC(-) modes. Table S-9: Detailed acquisition time in each MS mode. Figure S-1. Number of statistically significant features in urine analysis. (A) total feature number, (B) MS¹ annotated feature number, (C) MS² matched feature number. Figure S-2. Illustration of chromatograms in (A) Full-scan; (B) DDA; (C) DIA modes. Figure S-3. Number of total metabolic features as the function of parameter space (A) mass tolerance for constructing extracted ion chromatogram; (B) minimum peak height; (C) algorithm of chromatographic peak reconstruction.

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

Acknowledgment


This study was funded by University of British Columbia Start-up Grant (F18-03001) and Canadian Foundation for Innovation (CFI 38159). We sincerely appreciate the comments and suggestions received from MS experts of several major MS vendors, including Drs. Xuejun Peng (Bruker), John Vukovic (Waters), and Marcus Kim (Agilent). We also thank Ms. Alisa Hui for proof reading this manuscript.

Reference

- (1) Patti, G. J.; Yanes, O.; Siuzdak, G. *Nature reviews Molecular cell biology* **2012**, *13*, 263.
- (2) Wishart, D. S. *Nature reviews Drug discovery* **2016**, *15*, 473.
- (3) Warth, B.; Spangler, S.; Fang, M.; Johnson, C. H.; Forsberg, E. M.; Granados, A.; Martin, R. L.; Domingo-Almenara, X.; Huan, T.; Rinehart, D. *Analytical chemistry* **2017**, *89*, 11505-11513.
- (4) Beyer, B. A.; Fang, M.; Sadrian, B.; Montenegro-Burke, J. R.; Plaisted, W. C.; Kok, B. P.; Saez, E.; Kondo, T.; Siuzdak, G.; Lairson, L. L. *Nature chemical biology* **2018**, *14*, 22.
- (5) Rinschen, M. M.; Ivanisevic, J.; Giera, M.; Siuzdak, G. *Nature Reviews Molecular Cell Biology* **2019**, *1*.
- (6) Smith, C. A.; Want, E. J.; O'Maille, G.; Abagyan, R.; Siuzdak, G. *Analytical chemistry* **2006**, *78*, 779-787.
- (7) Pluskal, T.; Castillo, S.; Villar-Briones, A.; Orešič, M. *BMC bioinformatics* **2010**, *11*, 395.
- (8) Xia, J.; Sinelnikov, I. V.; Han, B.; Wishart, D. S. *Nucleic acids research* **2015**, *43*, W251-W257.
- (9) Benton, H. P.; Ivanisevic, J.; Mahieu, N. G.; Kurczy, M. E.; Johnson, C. H.; Franco, L.; Rinehart, D.; Valentine, E.; Gowda, H.; Ubhi, B. K. *Analytical chemistry* **2014**, *87*, 884-891.
- (10) Zhu, X.; Chen, Y.; Subramanian, R. *Analytical Chemistry* **2014**, *86*, 1202-1209.
- (11) Röst, H. L.; Rosenberger, G.; Navarro, P.; Gillet, L.; Miladinović, S. M.; Schubert, O. T.; Wolski, W.; Collins, B. C.; Malmström, J.; Malmström, L. *Nature biotechnology* **2014**, *32*, 219.
- (12) Siegel, D.; Meinema, A. C.; Permentier, H.; Hopfgartner, G. r.; Bischoff, R. *Analytical chemistry* **2014**, *86*, 5089-5100.
- (13) Schrimpe-Rutledge, A. C.; Codreanu, S. G.; Sherrod, S. D.; McLean, J. A. *Journal of the American Society for Mass Spectrometry* **2016**, *27*, 1897-1905.
- (14) Roemmelt, A. T.; Steuer, A. E.; Poetzsch, M.; Kraemer, T. *Analytical chemistry* **2014**, *86*, 11742-11749.
- (15) Nikolskiy, I.; Mahieu, N. G.; Chen, Y.-J.; Tautenhahn, R.; Patti, G. J. *Analytical chemistry* **2013**, *85*, 7713-7719.
- (16) Ivanisevic, J.; Zhu, Z.-J.; Plate, L.; Tautenhahn, R.; Chen, S.; O'Brien, P. J.; Johnson, C. H.; Marletta, M. A.; Patti, G. J.; Siuzdak, G. *Analytical chemistry* **2013**, *85*, 6876-6884.
- (17) Cajka, T.; Fiehn, O. *Analytical Chemistry* **2016**, *88*, 524-545.
- (18) Eliasson, M.; Rännar, S.; Madsen, R.; Donten, M. A.; Marsden-Edwards, E.; Moritz, T.; Shockcor, J. P.; Johansson, E.; Trygg, J. *Analytical Chemistry* **2012**, *84*, 6869-6876.
- (19) Libiseller, G.; Dvorzak, M.; Kleb, U.; Gander, E.; Eisenberg, T.; Madeo, F.; Neumann, S.; Trausinger, G.; Sinner, F.; Pieber, T.; Magnes, C. *BMC Bioinformatics* **2015**, *16*, 118.
- (20) Kale, N. S.; Haug, K.; Conesa, P.; Jayseelan, K.; Moreno, P.; Rocca-Serra, P.; Nainala, V. C.; Spicer, R. A.; Williams, M.; Li, X. *Current protocols in bioinformatics* **2016**, *53*, 14.13. 11-14.13. 18.
- (21) Zhou, J.; Li, Y.; Chen, X.; Zhong, L.; Yin, Y. *Talanta* **2017**, *164*, 128-136.
- (22) Wrona, M.; Mauriala, T.; Bateman, K. P.; Mortishire-Smith, R. J.; O'Connor, D. *Rapid Communications in Mass Spectrometry: An International Journal Devoted to the Rapid Dissemination of Up-to-the-Minute Research in Mass Spectrometry* **2005**, *19*, 2597-2602.
- (23) Hu, Y.; Cai, B.; Huan, T. *Analytical Chemistry* **2019**, *91*, 14433-14441.
- (24) Huan, T.; Li, L. *Analytical chemistry* **2015**, *87*, 7011-7016.
- (25) Shurubor, Y. I.; Paolucci, U.; Krasnikov, B. F.; Matson, W. R.; Kristal, B. S. *Metabolomics* **2005**, *1*, 75-85.
- (26) Rosenberger, G.; Koh, C. C.; Guo, T.; Röst, H. L.; Kouvonen, P.; Collins, B. C.; Heusel, M.; Liu, Y.; Caron, E.; Vichalkovski, A. *Scientific data* **2014**, *1*, 140031.
- (27) Hu, A.; Noble, W. S.; Wolf-Yadlin, A. *F1000Research* **2016**, *5*.
- (28) Doerr, A. *Nature methods* **2014**, *12*, 35.
- (29) Li, H.; Cai, Y.; Guo, Y.; Chen, F.; Zhu, Z.-J. *Analytical chemistry* **2016**, *88*, 8757-8764.
- (30) Yan, Z.; Yan, R. *Analytical chemistry* **2015**, *87*, 2861-2868.

- (31) Mullard, G.; Allwood, J. W.; Weber, R.; Brown, M.; Begley, P.; Hollywood, K. A.; Jones, M.; Unwin, R. D.; Bishop, P. N.; Cooper, G. J. *Metabolomics* **2015**, *11*, 1068-1080.
- (32) Calderón-Santiago, M. n.; Priego-Capote, F.; Luque de Castro, M. a. D. *Analytical chemistry* **2014**, *86*, 7558-7565.
- (33) Huan, T.; Forsberg, E. M.; Rinehart, D.; Johnson, C. H.; Ivanisevic, J.; Benton, H. P.; Fang, M.; Aisporna, A.; Hilmer, B.; Poole, F. L. *Nature methods* **2017**, *14*, 461.
- (34) Li, S.; Park, Y.; Duraisingham, S.; Strobel, F. H.; Khan, N.; Soltow, Q. A.; Jones, D. P.; Pulendran, B. *PLoS computational biology* **2013**, *9*, e1003123.
- (35) Pirhaji, L.; Milani, P.; Leidl, M.; Curran, T.; Avila-Pacheco, J.; Clish, C. B.; White, F. M.; Saghatelian, A.; Fraenkel, E. *Nature methods* **2016**, *13*, 770.
- (36) Wang, M.; Carver, J. J.; Phelan, V. V.; Sanchez, L. M.; Garg, N.; Peng, Y.; Nguyen, D. D.; Watrous, J.; Kapon, C. A.; Luzzatto-Knaan, T. *Nature biotechnology* **2016**, *34*, 828-837.
- (37) Yin, Y.; Wang, R.; Cai, Y.; Wang, Z.; Zhu, Z.-J. *Analytical Chemistry* **2019**, *91*, 11897-11904.
- (38) Samanipour, S.; Reid, M. J.; Bæk, K.; Thomas, K. V. *Environmental Science & Technology* **2018**, *52*, 4694-4701.
- (39) Nikolskiy, I.; Mahieu, N. G.; Chen, Y., Jr.; Tautenhahn, R.; Patti, G. J. *Analytical Chemistry* **2013**, *85*, 7713-7719.
- (40) Bern, M.; Finney, G.; Hoopmann, M. R.; Merrihew, G.; Toth, M. J.; MacCoss, M. J. *Analytical Chemistry* **2010**, *82*, 833-841.
- (41) Peckner, R.; Myers, S. A.; Jacome, A. S. V.; Egertson, J. D.; Abelin, J. G.; MacCoss, M. J.; Carr, S. A.; Jaffe, J. D. *Nature Methods* **2018**, *15*, 371-378.
- (42) Tsugawa, H.; Cajka, T.; Kind, T.; Ma, Y.; Higgins, B.; Ikeda, K.; Kanazawa, M.; VanderGheynst, J.; Fiehn, O.; Arita, M. *Nature methods* **2015**, *12*, 523.
- (43) Canterbury, J. D.; Merrihew, G. E.; MacCoss, M. J.; Goodlett, D. R.; Shaffer, S. A. *Journal of The American Society for Mass Spectrometry* **2014**, *25*, 2048-2059.
- (44) Myers, O. D.; Sumner, S. J.; Li, S.; Barnes, S.; Du, X. *Analytical chemistry* **2017**, *89*, 8689-8695.
- (45) Koopmans, F.; Ho, J. T.; Smit, A. B.; Li, K. W. *Proteomics* **2018**, *18*, 1700304.
- (46) Wu, I. L.; Turnipseed, S. B.; Storey, J. M.; Andersen, W. C.; Madson, M. *Rapid Communications in Mass Spectrometry* **2019**.
- (47) Eliasson, M.; Rännar, S.; Madsen, R.; Donten, M. A.; Marsden-Edwards, E.; Moritz, T.; Shockcor, J. P.; Johansson, E.; Trygg, J. *Analytical chemistry* **2012**, *84*, 6869-6876.

538 TOC graphical abstract



539

540

Figures




Figure 1. Schematic workflow for the comparison of full-scan, DDA and DIA modes in untargeted metabolomics.




Figure 2. Numbers of total metabolic features detected, MS¹ annotated and MS² matched metabolic features in urine samples analyzed by the three data acquisition methods in both (A) RP(+); (B) HILIC(-) analytical configurations.




Figure 3. RSD distribution of metabolic feature MS signal intensities in the analysis of (A) RP(+); (B) HILIC(-) analytical configurations.




Figure 4. Overlapped MS² features detected in DDA and DIA analysis in urine samples. (A) RP(+) analysis; (B) HILIC(-) analytical configurations.




Figure 5. Distribution of spectral similarity scores of DDA vs DIA in RP(+) and HILIC(-) analytical configurations.




Figure 6. Summary of the assessment of the three data acquisition methods.

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

569 **Tables**

570 **Table 1.** Number of confirmed metabolites in the standard mixture analyzed by three data
571 acquisition methods

Identified metabolites in standard mixture		Full-scan		DDA		DIA	
		No.	%*	No.	%	No.	%
StdMix(H)	RP(+)	66		63		66	
	HILIC(-)	125		120		125	
StdMix(L)	RP(+)	44	67	23	37	33	50
	HILIC(-)	57	46	43	36	46	37

572
573 *The percentage of detected metabolites in StdMix(L) compared to those in StdMix(H)